



Mageshape machine+Chair 2 in 1

# MAGESHAPE MIX



Muscle re-build of incontinent patients  
Body shaping and muscle building



# TIRED OF THE CONSTANT URGE TO GO?

Pelvic floor muscles are a layer of muscles that support the pelvic organs. Strong pelvic floor muscles provide control over the bladder and bowel. When these muscles are weakened, they do not fully support the internal organs, which may lead to difficulty controlling the release of urine.

## ARE YOU ?

Suffering from these afflictions  
Growing Need for Intimacy?



Dropped Wellbeing

Dissonance incouple's life



Coughing or Sneezing


Leaking urinewhile laughing




Decreased PelvicFloor  
Function  
Early AgingMenopause

# Traditional exercise

## KEGEL EXERCISE AND PELVIC EXERCISE



1. Preventing urinary incontinence and vaginal prolapse
2. Preventing uterine prolapse
3. Strengthening the pubococcygeal muscle



1. Male Prostate Pain
2. Benign prostatic hyperplasia
3. Prostatitis Magna Erectile Dysfunction

# New tech for improving urinary incontinence

## -Mageshape MIX

### What's Mageshape MIX?

A breakthrough treatment for incontinence and confidence. This unique technology revolutionizes women's intimate health and wellness by providing those suffering from incontinence with a completely non-invasive option.



# PRODUCT FEATURES



**Improve sexual sensitivity**  
Sensitivity and function



**Prevention and improvement**  
Stress urinary incontinence




**Postpartum recovery**  
Improve vaginal health



**Strengthen the pubic tail**  
Skeletal muscle (PC muscle)



**Increase blood in the pelvic floor**  
Flow and nerve function



**Increase the vagina after menopause**  
Wall thickness and lubricity

- ◆ Improve urine leakage and frequent urination.
- ◆ One course of treatment = 20,000 Kegel exercises
- ◆ Health care of the pelvic cavity
- ◆ Sex life is no longer affected
- ◆ No need to take off pants, expensive invasive treatment
- ◆ Both men and women can use



## How Does it Work?

Utilize high intensity focused electromagnetic techniques to stimulate the deep muscles of the pelvic floor and restore neuromuscular control. A single focused electromagnetic wave energy brings the greatest contraction on thousands of pelvic floor muscles, which is a very important muscle for incontinent patients. It makes the pelvic floor muscles perform hyper-contraction exercises (20,000 pelvic floor muscle contractions every 30 minutes).


-----




HIFEM stimulates deep layers with third-generation happiness ch




90% of people cannot exercise deep muscles properly



Unable to know whether the pelvic floor muscles are contracting correctly



The contraction frequency, intensity and amplitude are limited



100% contraction of deep pelvic floor muscles




Magnetic energy accurately focuses the pelvic floor muscle penetration depth up to 10 cm



Three-stage treatment mode is effective for training and strengthening



The strength and frequency of pelvic floor muscle contraction is limited



Ultra-limit contraction of pelvic floor muscles and maintain for a few seconds

# Why choose Mageshape Mix

## ➤ Non-invasive

Mageshape MIX treatment is non-invasive and does not require needles at all. HIFEM Energy stimulates the entire pelvic floor area and helps with restoring neuromuscular control.

## ➤ Improve quality of life


95% of patients have reported significant improvement in the quality of life.

## ➤ Comfortable

Mageshape MIX treatment is comfortable. You will experience a little tingling sensation throughout the 30 minutes at your pelvic area. No numbing cream is required.

## ➤ Effective

100% contraction of deep pelvic floor muscles, magnetic energy accurately focuses on the pelvic floor muscles with a penetration depth of 10 cm



30 = 20.000  
Mins Treatment Kegel exercise

30 MINS TREATMENT  
= 20.000 KEGEL EXERCISE





**2 treatment handles for  
targeted booy parts training**



Belly



Buttocks



Thighs



Belly




Legs



Arms

# HOW IT WORKS?



- » It uses HIEMT (High-Intensity Focused Electromagnetic Wave) technology to target motor nerve cells to stimulate muscle contraction with high intensity, which deeply reshapes the internal structure of muscles, strengthens the muscles, increases muscle density and volume.
- » Muscle extreme contraction requires a lot of energy supply, so the fat cells next to the muscle will be consumed, so it can help reduce fat effectively.
- » It makes body lines more prominent, "mermaid line" and "The vest line" clearer.
- » It can be said that "you can do exercise while lying down." In 30 minutes of treatment, it equals 20,000 workouts.

MINS TREATMENT

WORKOUTS

30 = 20.000

**16%** ↑ Average  
increase in  
muscle mass

**19%** ↓ Average  
fat reduction



# Operation System



# Operation System



# Operation System



# Android system

The image displays two screenshots of an Android application interface for a beauty salon management system.

**Left Screenshot (Registration/Profile Form):**

- Header: Profile picture of a woman.
- Form Fields:
  - Text input: ert
  - Text input: Phone/Email
  - Text input: age
  - Text input: height<sub>cm</sub>
  - Text input: weight<sub>kg</sub>
  - Gender selection: Male (♂) and Female (♀) icons.
- Action: CONFIRM button.
- Footer: ☐ "I have read and agreed the agreement"

**Right Screenshot (Dashboard/Calendar View):**

- Header: Back arrow, Date: 2024.10.11, User profile: ert.
- Summary Cards:
  - Login times: 2
  - Treatment counts: 0
- Date Selection:
  - YEAR: < 2024 > (2024.10.11)
  - MONTH: 06, 07, 08, 09, 10 (selected), 11, 12
- Calendar View (Table):

| SUN | MON | TUE | WED | THU | FRI | SAT |
|-----|-----|-----|-----|-----|-----|-----|
| 29  | 30  | 1   | 2   | 3   | 4   | 5   |
| 6   | 7   | 8   | 9   | 10  | 11  | 12  |
| 13  | 14  | 15  | 16  | 17  | 18  | 19  |
| 20  | 21  | 22  | 23  | 24  | 25  | 26  |
| 27  | 28  | 29  | 30  | 31  | 1   | 2   |
| 3   | 4   | 5   | 6   | 7   | 8   | 9   |

Beauticians administration • Patients administration

Clients information storage

Treatment parameter date memory system

Rental model

# Specification

|                        |                                       |
|------------------------|---------------------------------------|
| Magnetic wave(energy ) | 0-7tesla                              |
| Voltage                | 220V/110V                             |
| Out put power          | 3200W                                 |
| Frequency              | F1:1-10HZ    F2:1-100HZ               |
| Pulse width            | 300us                                 |
| Mode                   | Model-I( <b>Sport mode</b> )          |
|                        | Model-II( <b>Gentle mode</b> )        |
|                        | model-III( <b>Professional mode</b> ) |
| Screen                 | 10.4 inch                             |
| Machine size           | 830*715*780mm                         |
|                        |                                       |

# Frequently Asked Questions

**I have a low tolerance for pain! Is it painful?**

**A**

There's no pain at all. The experience is a tingling sensation with your pelvic floor muscles' contracting during the procedure. You may resume your daily activities immediately after the treatment.

**How long does the treatment take?**

**B**

A typical treatment would take 30-minutes. During the treatment, you will feel a slight tingling with no pain at all.

**How many sessions do I need?**

**C**

Always consult your provider for a tailored treatment plan. For optimal results, it's recommended to have a minimum 6 sessions, scheduled twice a week.

**How long will the result last?**

**D**

Our patients reported the results last up to 6 months. With regular kegel exercises, it can help to sustain the results past the 6-months mark!

**How about downtime and side effects?**

**E**

There is absolutely no downtime and side effects from the treatment. After the session, you can immediately return to your daily activities.